

# SHIELD<sup>me</sup>™

# FACE SHIELD



**SHIELD<sup>me</sup>™ FACE SHIELD** is used by general public, disinfection and sanitization team. Face Shield protects and prevents the entry of toxic liquid spatter or splash, harmful air borne microbes, and dust into your mouth, nose and eyes. It also prevents you from touching your eyes, mouth and nose subconsciously.

The lenses are made of PET and light transmittance material, transparent and clear enough to ensure not affecting the normal vision of the human eyes.

## Product Performance Features

- ✓ Overall ergonomic design, lightweight and more comfortable to wear.
- ✓ Made of polymer materials to effectively prevent foreign body impact and liquid spatter.
- ✓ The lenses are made of PET and light transmittance material with high clarity.
- ✓ Easy to adjust head band and suitable for all kinds of head size.



ANTI-IMPACT



ANTI-SPLATTER



LIGHTWEIGHT



LIGHT PENETRATION