

SHIELD^{me}™

NON-CONTACT IR THERMOMETER

SME-IRTM01

CE



SHIELD^{me}™ NON-CONTACT INFRARED THERMOMETER is a medical device for body / forehead temperature checking, super fast one-second reading with high measurement accuracy.

FEATURES:

With its innovative no-touch technology, this infrared thermometer is suitable for family members of all ages. It has smart built-in fever alarm system.

SPECIFICATION PARAMETER:

Product Name	Non-Contact IR Thermometer
Product Model	SME-IRTM01
Measurement range	32.0°C ~ 42.2°C (89.6°F - 108°F)
Accuracy	±0.2°C (35.0°C - 42.0°C) \ ±0.4°F (95.0°F - 107.6°F) ±0.3°C (±0.5°F) other range
Memory recall	32 Readings
Automatic shutdown	60s
Resolution	0.1°C
Operating conditions	16°C - 35.0°C (60.8°F to 95°F), 15% - 80% RH
Storage conditions	-20°C - +55°C (-4°F to 131°F), 15% - 93% RH

**We reserve the right to change specifications of our product without prior notice.*

info@shieldmeglobal.com

www.shieldmeglobal.com